IKA

designed for scientists

AOD 1.11 Control standard

5


/// DEUTSCH

9

/// ENGLISH

13

/// FRANÇAIS

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WORK INSTRUCTION FOR ION CHROMATOGRAPHIC AND TITRIMETRIC DETERMINATION OF IKA-CONTROL STANDARD AOD 1.11

EN

Product description:

IKA – Control standard AOD 1.11 for chlorine and sulphur

Shelf Life:

12 months after delivery if stored properly (check on the test certificate).

Requirements for proper storage:

- open the bottle only when taking the sample, close it again properly
- store below 10 °C
- protect from exposure to sunlight

Ion(s) sought: Chloride and Sulfate

Analysis by means of / with:

- 1) Autotitration / Titration processor, Titration of chloride with silver nitrate
- 2) Ion-chromatographic determination of chloride and sulfate by means of conductivity detection unit

The absorption solution is usually adapted to the eluent used with the Ion-Chromatograph

1 Liter stock solution is made up of:

21.2 g Na_2CO_3

6.30 g NaHCO_3

1 Liter absorption solution is made up of:

50 ml stock solution (see above)

25 ml 30% H_2O_2

7 - 8 NaOH – pellets

Calorific value of the sample:

Can be determined simultaneously when using C 7000, C 6000, C 5000, C 2000 or C 200 calorimeter

Combustion aid(s):

C 10 Acetobutyrate-capsules (IdNr.: 0000750000) may be used

O_2 - pressure: 30 bar

Description of the combustion residues:

Minimal residues of ash or soot residues in the crucible are possible.

Work instruction:

1. Sample preparation

From the homogeneous standard oil 0.2 - 0.5 g of the sample are filled directly into the quartz crucible or into an acetobutyrate-capsule (C 10).

The weighted-in quantity of the sample material is noted down for the later required calculation of the ion-content. The weight of the capsule is not important for the measurement as it contains neither halogens nor sulfur. The weight is important only during parallel calorific value determination, since the energy additionally introduced thereby must be deducted from the gross calorific value as foreign energy.

The capsule allows better sample handling and ensures that no volatile components of i.e. an unknown sample can escape. Furthermore, the capsule prevents splattering of the sample when igniting and supports the decomposition process through their additional energy provided.

The C 10 capsules may be used for the blank value determination as well.

2. Cleaning the decomposition vessel

It must be guaranteed that no contamination is present in the decomposition bomb, and the latter shall be cleaned in such a way that the parameters which may be examined are not affected.

Suitable cleaning agents are acetone, isopropyl alcohol and distilled water. The decomposition vessel is dried using a lint free cellulose cloth or the like and can be used for the decomposition.

We recommend performing blank value tests on a regular basis.

3. Preparing the decomposition vessel

First the cotton thread is bound onto the ignition wire. Then one adds 10 ml +/- 1ml of the absorption solution into the decomposition vessel (composition of the absorption solution see above).

Alternatively, it is possible to use 0,25 mol/l NaOH and 100 µl – 200 µl hydrogen peroxide (30%) as the absorption solution.

The quartz crucible is inserted into the crucible holder of the decomposition vessel. Then the cotton thread must be positioned inside the crucible in such a way, that it can ignite the acetobutyrate-capsule or oil directly.

The ignition insert is now placed into the decomposition vessel and screwed on tightly. Now the decomposition vessel is filled with 20 bar excess pressure and vented again two times by means of the filling station so that the entire nitrogen / oxygen of the room air (approximately 210 ml) is driven out. For the actual combustion the decomposition vessel is then filled with 30 bar oxygen (quality 3.5 or better) and is now ready for ignition. The gas-washing bottle of the venting station (see 4.a) is filled with 50 ml of the absorption solution.

4. The AOD 1 procedure

This section replaces point 4.a of this work instruction in the event that one uses the AOD 1 Decomposition System procedure.

A calorimeter for igniting the sample in the decomposition vessel is not needed for this procedure, instead one ignites the sample with an electrical impulse which is directly applied to

the electrodes of the decomposition vessel. One places the ignition adapter (AOD 1.13) onto the oxygen filled decomposition vessel (or places the vessel behind the safety wall AOD 1.3) and starts the combustion by pressing the igniter on the ignition box AOD 1.2. In any case, the ignition must be performed under safe conditions behind a safety wall and with enough safety distance. Please read the instruction manual of the AOD 1.1 decomposition vessel and AOD 1 System with care.

Successful ignition is characterized by the heating up of the decomposition vessel. After ignition the vessel can be cooled in a container filled with cold water. Once the vessel is cooled off, the incineration gases can be discharged carefully from the decomposition vessel via a venting station. The vessel may also be swiveled carefully to unite the absorption solution and condensed water on the wall of the decomposition vessel.

4.a Combustion in the calorimeter (C 7000 with C 7012 decomposition vessel)

One enters the required data for the measurement into the calorimeter, places the decomposition

vessel into the calorimeter and starts the measurement procedure. After the ignition the calorific value as well as the rise in temperature and calorific value are displayed. The decomposition vessel can be taken out of the calorimeter and be placed into the cooling unit of the system.

Once the vessel cooled off, the incineration gases are discharged in a slow, controlled manner from the decomposition vessel with the venting station C 7030 through a washing bottle filled with absorption solution.

Decomposition vessel	Suitable venting station	Calorimeter
C 5012	C 5030	C 5000, C 2000, C 200
C 6012	C 6030	C 6000
C 7012	C 7030	C 7000
AOD 1.1	C 7030	-

5. Rinsing the decomposition vessel for the potentiometric determination

The decomposition vessel is opened and thoroughly rinsed with distilled water. Care has to be taken, to ensure that all parts of the vessel are accounted for.

The decomposition solution as well as the absorption solution from the washing bottle are filled into a 250 ml flask and filled up to the calibration mark. The aliquot removed for the measurement (50.0 ml or 100.0 ml) is transferred into the titration cup.

The aliquot factor of 5 or 2.5 is noted and needed for calculation of the result. The titration cup is then filled approximately halfway with distilled water and 2 ml concentrated nitric acid is added.

6. Potentiometric measurement

The titration cup is then placed into the analysis processor and the corresponding procedure is started. After several minutes the result can be read off in milligrams. The halogen concentration can be calculated from this taking the originally weighted-in quantity and the aliquot factor into consideration.

7. Rinsing the decomposition vessel for analysis in the ion chromatograph

The procedure for this step is the same as described in paragraph 5., with the difference that

one uses the eluent of the ion-chromatograph instead of distilled water to rinse the decomposition vessel.

The calculation is performed by taking the originally weighted-in quantity into consideration as well as the flask contents in which the decomposition solution and the absorption solution were combined from the washing bottle.

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IKA India Private Limited
IKA England LTD.
IKA Vietnam Company Limited