2000005073

IKA®

C 48_012017

IKA® C 48

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Contents

	Page
Declaration of Conformity	7
Warranty	7
Warning symbols	7
Safety instructions	8
Correct use and description of the device	8
Unpacking	8
Installation	9
Commissioning	9
Maintenance and cleaning	10
Accessories	10
Technical data	10
List of spare parts-Illustration of spar parts	2

Declaration of conformity

We declare under our sole responsibility that this product corresponds to the regulations 2014/68/EU and conforms with the standards or standardized documents: EN ISO 12100.

Warranty

In accordance with **IKA**[®] warranty conditions, the warranty period is 12 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

Warning symbols

 DANGER	Indicates an (extremely) hazardous situation, which, if not avoided, will result in death, serious injury.
 WARNING	Indicates a potentially hazardous situation, which, if not avoided, can result in death, serious injury.
 CAUTION	Indicates a potentially hazardous situation, which, if not avoided, can result in injury.
 NOTICE	Indicates practices which, if not avoided, can result in equipment damage.

Safety instructions

For your protection

Read the operating instructions in full before starting up and follow the safety instructions.

- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.

Work with the device

DANGER

Before use, check the outlet pressure on the pressure reducer and compare it with the pressure displayed on the pressure gauge.

Oxygen as a compressed gas is oxidising; intensively aids combustion; can react violently with flammable materials.

- The filling time should not be set to less than 30 seconds. This will prevent combustible substances from being whirled around in the crucible.

- Insufficient ventilation can cause formation of explosive mixtures. Therefore, only use the device in well-ventilated areas.
- Check the seals for damage before each use.
- After finishing work, the oxygen supply on the pressure reducer must be closed.

For protection of the device

DANGER

Maintenance tasks should be performed without any excess pressure on the system.

The Valve needle (29) is held in place with safety varnish to keep it from being turned unintentionally. It should only be removed if there is no excess pressure on the system

- Keep lines and fittings through which oxygen flows free of grease.

Correct use and description of the device

Use

The **IKA**® C 48 oxygen filling station is only suitable for filling decomposition vessels C 7010, C 7012 and AOD 1.1 with oxygen.

Description of the device

- The 2 m pressure hose, included with delivery, can easily be laid to connect to the oxygen bottle.
- Therefore, the location of the C 48 filling station does not depend on the oxygen bottle being in the immediate vicinity.
- The filling pressure that is set and the pressure in the decomposition vessel can be read directly from the work station, thus preventing insufficient filling quantities because of not observing the minimum filling period.
- The filling time can be adjusted in increments.
- There is no need to be constantly occupied with the oxygen bottle because the automatic non-return valve only opens if a decomposition vessel is being filled.
- A carrier plate (screws on left-handed or right-handed) is used to prepare samples.

Unpacking

Unpacking

- Please unpack the device carefully.
- In the case of any damage a detailed report must be sent immediately (post, rail or forwarder).

Delivery scope

- C 48 oxygen filling station
- Holding plate
- O-rings
- Pressure hose (2 m)
- Operating instruction


Installation

Connecting the pressure hose

The pressure hose (43), included with delivery, is designed for a maximum pressure of 40 bar. It is 2 m long and easily be laid to connect to the oxygen bottle (smallest permissible bending radius 80 mm).


It should be connected at point "A (SW 10)".

Longer plastic lines are available on request.

There is an R1/4inch threading (SW 17) on the C 29 reducing valve (accessory) of the oxygen bottle to which the other end of the pressure hose can be connected.

An appropriate adapter is included for US pressure reducers.

Decomposition vessels

IKA® decomposition vessels C 1/8; C 7010, C 7011, C 7012 and AOD 1.1 are provided with an identification on the base or on the cap screw that contains information including the year of manufacturing number.

Assembly of the support plate

The support plate, included with delivery, can be screwed on right or left. The cover of the decomposition vessel can be set down on the support plate. This makes it easier to load decomposition vessels.

Commissioning


The required filling pressure -usually 30 bar -can be set on the pressure reducer valve of the oxygen bottle. The adjusted pressure can also be read on the pressure gauge (11).

The automatic non-return valve (33) integrated into the oxygen filling station prevents oxygen from flowing out when this is not intended.


When the task is complete, the main valve on the oxygen supply should be closed.

The decomposition vessel should be placed on the holding plate (38 or 50). By loosening the screw (3), the holding plate (38 or 50) can be adjusted so that the filling head (32) fits exactly onto the oxygen line valve of the decomposition vessel.

All you need to do then to fill up the decomposition vessel is to swing the hand lever downwards. The increase in pressure can be observed on the pressure gauge.

The factory-set filling time can be changed by turning on the valve needle (29). Turning the safety valve counter-clockwise results in a longer filling time.

The filling time should not be less than 30 seconds. This will prevent combustible substances from being whiled around in the crucible.

The valve needle (29) is held in place with safety varnish to keep it from being turned unintentionally. It should only be turned if there is no excess pressure on the system.

Maintenance and cleaning

The C 48 oxygen filling station requires no maintenance to operate. It is subject only to the natural aging of the components and their statistical failure rate.

When ordering spare parts, please indicate the manufacturing number on the rating plate, the device type and the item number and description of the spare parts.

Please send in only units that have been cleaned and that are free of hazardous substances for repairs. Use only cleaning agents approved by **IKA®** to clean **IKA®** devices. To clean the residues below, use the corresponding cleaning agents:

Dye, coloring	Isopropanol
Construction materials	Water containing detergent
Cosmetics	Water containing detergent/isopropanol
Foodstuffs	Water containing detergent
Combustible materials	Water containing detergent


For materials that are not listed above, please inquire with **IKA®**. Wear protective gloves to clean the units. Electrical devices must not be placed in the cleaning agent to clean them. Before using another cleaning or decontamination method than the one recommended by the manufacturer, the user must make certain by consulting with the manufacturer that the intended method will not damage or destroy the unit.

It should be expected that frequent use will result in wear on the O rings in the filling head.

The filling head is screwed in place and can be unscrewed after the pressure hose is removed. The lower O ring (21), which is subject to particularly high wear and tear, can be reached without disassembling the pressure hose after the cap screw (34) has been removed.

Item 23 Seal for the non-return valve

Item 21 Seal for the decomposition vessel

Repair

In case of repair the device has to be cleaned and free from any materials which may constitute a health hazard.

For this you should use the downloaded print of the „**Certificate of compliance**“ form from the **IKA®** website www.ika.com.

If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

Accessories

- C 29 Pressure gauge, oxygen
- Holding plate: for inserting decomposition vessels C 7010, C 7011, C 7012, C 7010 Fz C 7012 FZ and AOD 1.1

Technical data

Permissible operating pressure	bar	max. 40
Pressure gauge	bar	0-50
Typical filling pressure		30 steplessly adjustable
Pressure hose	m	2
Filling speed	sec	> 30 (preset to 30 bar)
Current filling pressure		readable
Oxygen valve		automatic non-return valve
Dimensions (B x T x H)	mm	200 x 260 x 450
Weight	kg	6

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IKA®-Werke GmbH & Co.KG