

## Immunofluorescence Staining Using the $\mu$ -Slide 18 Well

This Application Note presents a simple protocol for the cultivation, fixation, and staining of adherent cells using ibidi's  $\mu$ -Slide 18 Well. In this example, we cultivated human endothelial cells, fixed them with paraformaldehyde, and stained the F-actin cytoskeleton and alpha-tubulin expressing microtubules. The nuclei were counterstained with DAPI.

### Related Topics:

[Application Note 09 Immunofluorescence Staining with  \$\mu\$ -Slide VI<sup>0.4</sup>](#)

[Application Note 16 Immunofluorescence Staining with  \$\mu\$ -Slide 8 Well](#)

### Keywords:


Immunofluorescence, Fixation, Staining, Mounting, Microscopy, Cell Culture, HUVECs, F-actin, alpha-tubulin, DAPI, Chambered coverslip

### Table of contents

0. General Information.....	1
1. Cultivation.....	2
2. Fixation, Permeabilization & Blocking.....	2
3. Staining.....	2
4. Imaging.....	3
5. Results.....	3

### 0. General Information

The protocol consists of four main steps:



In this application note, the following material is used:

Material	Manufacturer	Catalog number
μ-Slide 18 Well ibiTreat	ibidi	81816
HUVEC	various	various
4% Paraformaldehyde	Sigma-Aldrich	HT5011
0.1% Triton® X-100 (diluted in PBS)	Alfa Aesar	A16046
1% Bovine Serum Albumin (diluted in PBS)	Sigma-Aldrich	A1470
LifeAct-TagGFP2 Protein	ibidi	60112
Monoclonal anti-alpha-Tubulin antibody, mouse	Sigma-Aldrich	T5168
Anti-mouse IgG-Atto594	Sigma-Aldrich	76085
ibidi Mounting Medium with DAPI	ibidi	50011
Optional: ibidi Immersion Oil	ibidi	50101
Fluorescence microscope (inverted) with appropriate filter sets		

## 1. Cultivation

- Unpack the ibidi μ-Slide 18 Well ibiTreat (#81816) under sterile conditions and put it on a μ-Slide Rack (#80003).
- Prepare the cell suspension ( $1 \times 10^5$  cells/ml) and apply 100 μl into each well of the μ-Slide 18 Well.
- Cover the chambers with the supplied lid.
- Cultivate overnight in a humid cell culture incubator (37°C, 5% CO<sub>2</sub>). For longer cell cultivation, we recommend a medium exchange after a few days.

## 2. Fixation, Permeabilization & Blocking

- Aspirate the cell culture medium from the wells using a cell culture aspiration device.
- Wash cells with Dulbecco's PBS by slowly applying 100 μl into each well.
- Fix cells with ~100 μl of 4% paraformaldehyde for 20 min.
- Wash cells twice with PBS by slowly applying 100 μl into each well.
- Remove the liquid and apply ~100 μl of 0.1% Triton® X-100 in PBS. Incubate for 10 min.
- Wash cells with PBS by slowly applying 100 μl into each well.
- Remove the liquid and apply 100 μl 1% BSA blocking solution in PBS. Incubate for 20 minutes.
- Wash cells with PBS by slowly applying 100 μl into each well.

## 3. Staining


- Prepare your staining and antibody solutions.
- Remove all liquid from the wells using a cell culture aspiration device.
- Apply 100 μl of primary antibody (anti-Tubulin 1:1000) diluted in PBS into each well and incubate overnight at 4 °C.

- Wash the cells once with PBS for 10 minutes.
- Apply 100  $\mu$ l of secondary antibody mix (Anti-mouse IgG-Atto594 1:500 and LifeAct-TagGFP2 Protein 30  $\mu$ g/ml) diluted in PBS into each well and incubate for 3 hours at room temperature in the dark.
- Wash the cells twice with PBS for 10 minutes.
- Aspirate the PBS and add  $\sim$ 100  $\mu$ l ibidi Mounting Medium with DAPI into each well. Use the dropper bottle to add the Mounting Medium dropwise.

#### 4. Imaging

- Observe the cells under a fluorescence microscope with appropriate filter sets and optionally with ibidi Immersion Oil (#50101).
- Optionally, overlay images to create a merged image.

#### 5. Results



*HUVEC (Human umbilical vein endothelial cell) 24 hours after seeding in  $\mu$ -Slide 18 Well, ibiTreat, objective lens 60x, oil immersion.*

*Green: F-Actin (LifeAct-TagGFP2 protein)*

*Red: Tubulin (monoclonal anti-alpha-Tubulin antibody, mouse + Anti-mouse IgG-Atto594)*

*Blue: Nuclei (DAPI staining using ibidi Mounting Medium with DAPI)*