



---

PRODUCT INFORMATION

## Rapid Blood Genomic DNA Extraction Kit

*Product information for BT4782:*

### Kit Contents

Components	BT4782, 50 Preps
Buffer TBP	40 ml
Universal Buffer Digestion	10 ml
Buffer PR	4 ml
Proteinase K	1.2 ml
TE Buffer	50 ml
Protocol	1

### Storage

Transportation at ambient temperature, Store at 4°C, Valid for 1 year.

### Introduction

The kit allows simple and fast isolation of high quality genomic DNA from fresh and anti-coagulated blood. The entire procedure takes approx. 15 minutes only. Up to 5 ml of blood sample can be treated for each mini-preps. Average DNA yields are 20–60 µg per ml of the whole blood sample. Purified DNA can be used in a wide range of downstream applications.



---

## Features

- ✓ High quality.  $OD_{260}/OD_{280}$  of purified DNA is between 1.8~2.0, and can be directly used for downstream applications, such as PCR, Restriction Enzyme Digestion, Southern blotting.
- ✓ Simple procedure. Take 15 min to complete the whole procedure.
- ✓ High Yield. The yield of DNA is generally at 20~60  $\mu$ g of 1 ml whole blood, and up to 5 ml of whole blood can be used in a single prep.
- ✓ Easy to scale up.
- ✓ Non toxic.

## Procedures

1. Transfer 300 $\mu$ l anti-coagulated blood to a clean 1.5 ml tube. Add 600 $\mu$ l Buffer TBP, mix well and incubate for 1 min at room temperature.

**Note 1:** For frozen blood sample, pre-warm sample at 37°C for 3 minutes.

**Note 2:** For a fowl blood sample, transfer 100 $\mu$ l of the blood sample to 1.5 ml microtube, go to step 4.

2. Centrifuge at 12,000 x g for 1 min at room temperature, discard supernatant carefully.

3. Add 500 $\mu$ l TE Buffer, suspend the deposit gently but thoroughly, Centrifuge at 12,000 x g for 1 min at room temperature, discard supernatant carefully.

**Note:** Repeat this step once if necessary.

4. Add 180 $\mu$ l Buffer Digestion and 20 $\mu$ l Proteinase K, mix thoroughly. Incubate at 56°C for 20~30 min.

**Note:** To obtain RNA-free DNA, add 20 $\mu$ l RNase A solution (20 mg/ml) to the tube, mix thoroughly and incubate at room temperature for 5 minutes.



- 
5. Add 60 $\mu$ l of Buffer PR, mix by inverting 5~10 times, incubate at -20°C for 20 minutes.
  6. Centrifuge at 12,000 x *g* for 5 minutes at room temperature. Transfer the supernatant to a new 1.5 ml tube.
  7. Add equal volume of isopropanol (approx 0.15~0.25 ml) to the solution, mix well by inverting 5 times. Incubate at room temperature for 2~5 minutes. Centrifuge at 12,000 x *g* for 5 minutes, discard the supernatant carefully.
  8. Add 1 ml of pre-cooled 75% ethanol to the tube, mix well by inverting 10 times. Centrifuge at 12,000 x *g* for 1 minute, discard the supernatant.
  9. Repeat the Step 8 once.
  10. Air-dry the pellet at room temperature with the lid open for 2~5 minutes.

**Note:** Don't over dry.

11. Add 50~200 $\mu$ l of TE buffer to dissolve DNA pellet. Keep at 4°C for a couple hours until DNA pellet is completely dissolved. Purified DNA is ready for use. Or store at -20°C for long term storage.



**Bio Basic Inc.**

Quality-Affordable Research

---

**PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC  
RESEARCH ONLY!  
NOT INTENDED FOR HUMAN OR ANIMAL USE!**

Please visit [www.biobasic.com](http://www.biobasic.com)



**Quality-Affordable Research**