
20000006429

C 5030_7030_102015

C 5030

C 7030

Betriebsanleitung	DE	3
Ursprungssprache		
Operating instructions	EN	8
Mode d'emploi	FR	13
Руководство пользователя	RU	18

IKA®-Werke, Germany
Reg. No. 004343

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EN

Contents

	Page		Page
Declaration of conformity	2	Setting up of the device	9
Contents	8	Function	10
Explanations of symbols	8	Maintenance and cleaning	11
Safety instructions	8	Fault table	12
Unpacking	9	Warranty	12

Explanations of symbols

General hazard

This symbol identifies information **which is of absolute importance for the safety of your health**. Non-compliance may lead to impairments to health and injury.

This symbol identifies information which is important for proper operation and handling of the decomposition vessel. Non-compliance may result in imprecise measurements or damage to the venting station.

Safety instructions

Read the operating instructions in full before starting up and follow the safety instructions.

Keep the operating instructions in a place where they can be accessed by everyone.

Ensure that only trained staff work with the appliance.

Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.

Furthermore, toxic combustion residue in the form of gases, ashes or condensation, for example, is possible in the inner wall of the washing bottle.

Observe the accident prevention requirements applicable to the activity and the work station.

Wear your personal protective equipment.

When handling combustion samples, combustion residue and auxiliary materials, please observe the relevant safety regulations. The following materials, for example, could pose a risk:

- corrosive,
- highly flammable,
- explosive,
- bacteriologically contaminated,
- toxic.

Combustion gases are hazardous to health, therefore the venting hose must be connected to a suitable gas cleaning system or extraction system.

The C 5030 venting station is used for controlled venting of **IKA**[®] decomposition vessels C 5010, C 5011 and C 5012 in conjunction with the **IKA**[®] C 5000 control / duo control calorimeter system in the "decomposition" operating mode.

The C 7030 venting station is used for controlled venting of **IKA**[®] decomposition vessels C 7010, C 7011, C 7012 and AOD 1.1.

The venting speed may be infinitely adjusted via a throttle valve by using a manual adjusting screw. Please note that the throttle valve is **closed** at the beginning of the venting operation.

Do not operate without the ground clamp included in the delivery quantity!

The maximum permitted operating pressure is 40 bar (at the vent header).

Please observe the general laboratory guidelines when handling glass vessels..


Unpacking

• Unpacking

- Please unpack the device carefully
- In the case of any damage a report must be sent immediately (post, rail or forwarder).


• Delivery quantity

- 1 Gas washing bottle holder
 - 2 Gas washing bottle
 - 3 GL screw joint
 - 4 KC 29 ground clamp
 - 5 Venting head
 - 6 Drainage hose
 - 7 Pressure tubing
- A/F 10 fork spanner
Set of O-rings
Operating instructions



• Delivery quantity


- 1 Base plate, erect.
 - 2 Gas washing bottle
 - 3 GL screw joint
 - 4 Pressure tubing
 - 5 KC 29 ground clamp
 - 6 Venting grip C 7030
 - 11 Drainage hose
- A/F 10 single end spanner
A/F 8/10 twin end spanner
O-Ring 5x2
O-Ring 8x2,5
O-Ring 5x1,5
Operating instructions

Setting up of the device

Place the venting station on a horizontal surface and insert the gas washing bottle into the gas washing bottle holder until it comes to a stop. Connect the drainage hose to the throttle valve and the GL screw joint of the washing bottle as shown in the diagram. The tubing is fastened by closing the GL screw joint.

C 5030

Function

Pour an absorption solution that is suited to your purposes into the gas washing bottle and **secure the ground joints by means of the ground clamp**. Close the throttle valve by slightly turning the adjusting screw in a clockwise direction.


Note:

The throttle valve does not function as a shut-off valve. There may be some formation bubbles in the system when the adjusting screw is hand-tightened with maximum operating pressure in the decomposition vessel.

WARNING

Now fasten the pressure tubing to the venting head and the throttle valve. All of the screw joints are tightened by means of the A/F 10 fork spanner which is included in the delivery quantity.

C 7030

C 5030

Remove the decomposition vessel from the calorimeter and place it into the venting station.

Push the venting head on to the decomposition vessel until it comes to a stop. The valve of the decomposition vessel is opened by pressing down the self-locking eccentric lever.

The venting head cannot be separated from the decomposition vessel in this state.

Check to make certain that the ground clamp fits correctly onto the gas washing bottle and slowly open the throttle valve.

Take note of any bubble formation and correct in order to achieve the venting speed desired.

The venting head is removed when the decomposition vessel is completely expanded (lever in upward position) and the throttle screw is closed again.

WARNING

The gas washing bottle may be damaged if venting is performed too quickly.


C 7030

Remove the decomposition vessel from the calorimeter/cooler or remove the ignition device and place it on the venting station.

The venting screw on the handle should be screwed out to about 6 or 7 mm. Place the venting handle on the decomposition vessel and lock it in place by turning in a clockwise direction. As the screw lowers (being turned to the right) the lock valve of the decomposition vessel is opened.

In this state, the venting handle cannot be separated from the decomposition vessel.

Check to make certain the ground clamp is fitted snugly in the gas washing bottle and slowly open the throttle valve. Observe the formation of bubbles and correct until the desired venting speed is achieved. When the pressure has been entirely released from the decomposition vessel, unscrew the lock by rotating the entire handle to the left. Close the throttle screw on the venting valve.

⚠ WARNING The gas washing bottle may be damaged if venting is performed too quickly.

Maintenance and cleaning

C 5030

Clean the capillary boring of the venting flask from time to time through expansion with a decomposition vessel filled with pure oxygen-filled without a gas washing bottle attached.

Check the condition of the sealing washer on the filling flask and replace it if necessary.

Flask replacement / removing the pressure spring:

The venting head is disassembled by loosening the three countersunk screws (Item 80) located on it. The flask can then be completely removed.

C 7030

Clean the venting handle from time to time by releasing the gas pressure in a decomposition vessel filled with pure oxygen without any washing bottle attached.

Check the condition of the O-ring seals on the throttle screw of the valve as well as on the venting handle. These seals are subject to natural wear and tear. If the system does not seal completely and the seals show signs of mechanical wear, they must be replaced.

Only use cleansing agents which have been recommended by IKA®.

Use to remove:

Dyes	Isopropyl alcohol
Construction materials	Water containing surfactant, Isopropyl alcohol
Cosmetics	Water containing surfactant, Isopropyl alcohol
Foodstuffs	Water containing surfactant
Fuels	Water containing surfactant

- Do not allow moisture to get into the appliance when cleaning.
- Wear protective gloves when cleaning the devices.
- Please consult with IKA® before using any cleaning or decontamination methods, other than those recommended here.

Ordering spare parts

When ordering spare parts, please give:

- Device type
- serial number, see type plate
- Item number and designation of the spare part, see www.ika.com.

Repair

Please send instrument in for repair only after it has been cleaned and is free from any materials which may constitute a health hazard.

For this you should request the “**Decontamination Certificate**” from **IKA®**, or use the download printout of it from the **IKA®** website **www.ika.com**.

Return the instrument in its original packaging. Storage packaging is not sufficient. Also, please use suitable shipping package materials.

Fault table

C 5030

Error: Decomposition valve is not expanded.
Cause: Flask capillary is deformed.
Remedy: Replace flask.

Error: Leakage on venting head when expanding.
Cause: Sealing washer missing or defective.
Remedy: Replace sealing.

Error: Flask in venting head does not retract.
Cause: Flask boring is dirty or the pressure spring is defective.
Remedy: Clean or replace spring.

Error: Leakage on throttle valve
Cause: O-ring on throttle screw defective.
Remedy: Replace O-ring.

C 7030

Error: Leakage on venting head when expanding.
Cause: Sealing washer missing or defective.
Remedy: Turn the venting screw counterclockwise until it comes completely out and replace the O-ring (Pos. 23); or pry O-ring (Pos. 22) out of the venting handle and set a new ring in place.

Error: Leakage on the throttle valve
Cause: O-ring on throttle screw defective.
Remedy: Turn the throttle screw counterclockwise until it comes completely out and replace the o-ring (Pos. 45).

Error: Throttle valve exhibits high leakage rate when it is closed.
Cause: Needle valve of the throttle screw is excessively ty or corroded.
Remedy: Screw throttle screw completely out. Clean or replace.

Warranty

In accordance with **IKA®** warranty conditions, the warranty period is 12 months. For claims under the warranty please contact your local dealer. You may also send the machine directly to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

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IKA® - Werke GmbH & Co.KG