

## Generation and Dynamic Culture of L929 Spheroids in the $\mu$ -Slide Spheroid Perfusion

This Application Note is an example protocol for creating multicellular spheroids in the  $\mu$ -Slide Spheroid Perfusion with subsequent flow application. After formation of spheroids from a suspension of the murine fibroblast cell line L929, perfusion is applied with the ibidi Pump System. This ensures the optimal nutrition of spheroids during long-term cultivation.

### Related Documents:

- [Instructions  \$\mu\$ -Slide Spheroid Perfusion](#)
- [Instructions ibidi Pump System](#)

### Keywords:

Spheroids, organoids, 3D aggregates, L929 cells, long-term culture, Bioinert, passivation, flow, perfusion, pump, microscopy

### Material:

- $\mu$ -Slide Spheroid Perfusion Bioinert (80350, ibidi, Germany)
- Murine fibroblast cell line L929 (ACC 2, DSMZ, Germany)
- Cell culture medium (RPMI-1640 (21875034) with FCS (10270106), Gibco)
- Accutase (A1110501, Gibco)
- ibidi Pump System (10902, ibidi, Germany)
- Perfusion Set BLUE, 15 cm, ID 0.8 mm (10961, ibidi, Germany)
- Standard cell culture equipment (sterile working bench, cell culture incubator, culture flasks, PBS, etc.)

**Important Note:** Equilibrate all required materials, such as  $\mu$ -Slides, culture medium, and tubing (Perfusion Sets), **overnight** inside the incubator at 37°C and 5% CO<sub>2</sub>. This is essential for keeping air bubbles from emerging over time.

## 1. Cell Preparation & Seeding

- Cultivate L929 cells in culture medium.
- Treat the cells with Accutase for 1–2 minutes for detachment.
- Harvest the cell suspension.
- Centrifuge the cell suspension and dilute it in culture medium; the amount depends on the required cell concentration.
- Count the cells and adjust to a concentration of  $5 \times 10^5$  cells/ml.
- Prepare the  $\mu$ -Slide according to the instructions.
- Seed the cells into the channel of the  $\mu$ -Slide. Do not forget to seed cells for a static control.
- Incubate over night at 37°C and 5% CO<sub>2</sub> for spheroid formation.

## 2. Spheroid Formation


- Control the spheroid formation under the phase contrast microscope.
- Optional: after the spheroids have formed, wash 1x with fresh culture medium.
- Incubate for one hour at 37°C and 5% CO<sub>2</sub> before starting the perfusion.

## 3. Perfusion Experiment


- Prepare the  $\mu$ -Slide, the ibidi Pump System, and the Perfusion Set for the flow connection according to the instructions.
- To remove air bubbles from the system, let it run 1–2 hours before connecting the  $\mu$ -Slide.
- Connect the  $\mu$ -Slide to the tubing and the pump system.
- Start the perfusion experiment with 5 mbar resulting in a flow rate of 0.75 ml/min.
- Determine the flow rate. If necessary, adjust the pressure to create the desired flow rate.
- For the static control, perform a medium exchange every two days according to the instructions of the  $\mu$ -Slide Spheroid Perfusion.

## 4. Results

Spheroid formation starts directly after seeding. During long-term culture, the application of flow leads to stronger proliferation of spheroids, resulting in an increase of the spheroid diameter due to optimal nutrition.



*L929 fibroblasts show spheroid formation in the  $\mu$ -Slide Spheroid Perfusion, Bioinert, days 1–14, seeding concentration  $5 \times 10^5$  single cells/ml. Top: perfusion with the ibidi Pump System, 0.75 ml/min. Bottom: no perfusion, medium exchange every second day. Phase contrast microscopy, 10x objective lens, well diameter 800  $\mu$ m.*



*Comparison of spheroid formation of L929 fibroblasts in the  $\mu$ -Slide Spheroid Perfusion, Bioinert, seeding concentration  $5 \times 10^5$  single cells/ml. Static: no perfusion, medium exchange every second day. Perfusion: perfusion with the ibidi Pump System, 0.75 ml/min.*