



PRODUCT INFORMATION

Rapid Bacteria RNA Isolation Kit

Product information for BS8625:

Kit Contents

Components	BS8625, 50 Preps
Buffer Rlysis-B	50 ml
RNase-free Water	5 ml
Protocol	1

Storage

Store at room temperature. The kit is valid for 1 year at 4°C.

Introduction

This kit is designed for preparation of high quality total RNA from bacteria cells. 20 µg total RNA can be purified from 5 x 10⁷ bacteria cells using this kit. Purified RNA is ready for most downstream applications such as RT-PCR, Northern Blotting, Poly (A) purification, nuclease protection and in vitro translation.

NOTE: Care must be taken when working with RNA. It is important to maintain an RNase-free environment starting with RNA sample preparation and continue through purification and analysis. Use RNase free tubes, tips, gels. Wear gloves at all times.



Features

- ü Fast. Using fast lysis buffer, the whole procedure takes less than 40 minutes.
- ü High Quality of RNA. Purified RNA has an OD₂₆₀/OD₂₈₀ ratio of 1.9-2.0.
- ü Easy to scale up

Materials Supplied by User:

Microcentrifuge capable of at least 12,000 × g
RNase-Free pipets and pipet tips
Vortexer
RNase-Free Ethanol (96-100%)
RNase-Free Microcentrifuge tubes (1.5 ml or 2 ml)

Procedures

1. Sample Preparation.
 - A. Gram-negative bacteria (*E. coli*, streptococcal, pneumococcal, etc.)
 - a. Transfer logarithmic phase culture (about 2 x 10⁹ cells) into centrifuge tube and centrifuge at 10000 x g for 30 seconds, discard supernatant.
 - b. Add 100 µl lysozyme solution (400 µg/ml lysozyme in RNase-free Water. NOT supplied in the kit), suspend thoroughly and incubate at 37°C for 5 minutes.
 - B. Gram-positive bacterial (golden staphylococcal, orynebacteriadiphtheriae, etc.)
 - a. Transfer logarithmic phase culture (about 2 x 10⁹ cells) into centrifuge tube and centrifuge at 10000 x g for 30 seconds, discard supernatant.



Bio Basic Inc.

A world leader in serving science

- b. Add 100 μ l lysozyme solution (3 mg/ml lysozyme in RNase-free Water. NOT supplied in the kit), suspend thoroughly and incubate at 37°C for 5 minutes.
2. Using RNase-free pipet tips, add 1 ml Buffer Rlysis-B and mix by inverting immediately.
3. Incubate at room temperature for 5 minutes to make sure the cells are completely lysed.
4. Add 200 μ l chloroform to the tube, mix by inverting.
5. Centrifuge at 12,000 \times g for 5 minutes at 4°C. Transfer the supernatant to a new RNase-free 1.5 ml tube.
6. Add 1/3 volume of ethanol to the tube, vortex for 30 seconds.

Note: Incubating at -20°C for 5-10 minutes can improve the RNA yield.
7. Centrifuge at 12,000 \times g for 5 minutes at 4°C, discard the supernatant carefully.
8. Add 1 ml of RNase-free 75% ethanol to the tube, mix well by inverting 10 times. Centrifuge at 12,000 \times g for 1 minute, discard the supernatant.
9. Repeat the Step 8 once.
10. Air-dry the pellet at room temperature with the lid open for 2~5 minutes.

Note1: This step is very important, residual ethanol in RNA will interfere some downstream applications.

Note2: Don't over dry.
11. Add 30~50 μ l of RNase-free Water to dissolve RNA pellet. The purified RNA is ready for use. Or keep at -70°C for long term storage.



Bio Basic Inc.

A world leader in serving science

**PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC
RESEARCH ONLY!
NOT INTENDED FOR HUMAN OR ANIMAL USE!**

Please visit www.biobasic.com



A world Leader in Serving Science