

# IKA

designed for scientists

EN



MORE THAN JUST DISPERSING | IKA ULTRA-TURRAX®  
Tube Drive Power control

# ULTRA-TURRAX® Tube Drive P control

/// 4 in 1 without any cross-contamination

The new IKA UTDD power control disperser gets its name from an even more powerful motor. The disperser assures the best stirring, dispersing, grinding and homogenisation results and has been devised in particular for the safe processing of infectious, toxic or odour-intensive samples.



## Disposable tubes

The sealable disposable tubes have been optimised for your application. They reduce the amount of cleaning required and rule out any risk of cross-contamination.

## Programmable test conditions

Time, energy, volume: You specify and systematically program the test conditions and store frequently used parameters in your library.

## Adjustable change in direction of rotation

When, and how often, the direction of rotation should be changed is very easily set using the clear and multilingual menu.

## Turbo button

Even more intensive mixing, dispersing or grinding at your fingertips: All of this is achieved with the turbo button. Over short periods, rotational speeds of up to 8,000 rpm can be achieved.



ULTRA-TURRAX® Tube Drive P control  
Ident. No. 0025005981