

EPPENDORF pipettes		Reference 2 single channel						Reference 2 single channel			Research Plus single channel										Research Plus multichannel						
HTL tips	Tip volume [µl]	Cat number		0.1-2,5	0.5-10	2-20	10-100	30-300	100-1000	0.5-10	10-100	30-300	0.1-2,5	0.5-10	2-20	5-100	10-100	20-200	30-300	100-1000	500-5000	0.5-10	10-100	30-300			
NON FILTER	0.5-10	301-03-052, 301-22-052		•						•			•									•					
		301-03-102, 301-22-102		•										•									•				
		301-03-152, 301-22-152		•											•									•			
		301-03-351		•											•									•			
		301-79-102		•											•									•			
	1-200	301-02-052, 301-06-052				•																					
		301-02-152				•																					
		301-02-202, 301-06-352				•																					
		301-06-252				•																					
		301-79-052				•																					
1-300	301-16-052				•																						
	301-16-102				•																						
100-1000	301-01-202								•																		
	301-01-302									•																	
500-5000	301-01-352								•																		
	301-71-052																					•					
FILTER	0.5-10	302-05-102, 302-26-152		•																							
	0.5-20	302-03-102				•																					
	1-100	302-08-102					•																				
	1-200	302-04-102																									
	1-300	302-12-102																									
	100-1000	302-01-102																									

• compatible

• incompatible

• compatible, but the ejector does not eject the tip

November 2020