

GILSON pipettes		Pipetman Classic single channel					Pipetman G single channel					Pipetman L single channel					Pipetman L multichannel				Pipetman Neo single channel					Pipetman Neo multichannel							
HTL tips	Tip volume [µl]	Cat number	0.2	1-10	2-20	20-100	50-200	200-1000	0.2	1-10	2-20	10-100	20-200	100-1000	0.5-10	2-20	20-200	20-300	0.2	1-10	2-20	10-100	20-200	100-1000	0.2	1-10	2-20	10-100	20-200	100-1000	2-20	20-200	
	0.5-10	301-03-052, 301-22-052	•						•	•										•	•					•	•						
	0.5-10	301-03-102, 301-22-102	•	•					•	•										•	•					•	•						
	0.5-10	301-03-152, 301-22-152	•	•					•	•										•	•					•	•						
		301-03-351	•	•					•	•										•	•					•	•						
		301-79-102	•						•											•						•							
		301-02-052, 301-06-052		•						•						•																•	
		301-02-152		•						•						•																•	
	1-200	301-02-202, 301-06-352		•						•						•																•	
		301-06-252		•						•						•																•	
		301-79-052		•						•						•																•	
	1-300	301-16-052		•						•						•																•	•
		301-16-102		•						•						•																•	•
	100-1000	301-01-202						•						•																		•	•
		301-01-302						•						•																		•	•
		301-01-352						•						•																		•	•
	0.5-10	302-05-102, 302-26-152	•																														
	0.5-20	302-03-102		•						•						•																•	
	1-100	302-08-102		•						•						•																•	
	1-200	302-04-102		•						•						•																•	•
	1-300	302-12-102		•						•						•																•	•
	100-1000	302-01-102						•						•																		•	•

• compatible

• incompatible

November 2020