

SOCOREX pipettes		Accura Manual single channel								Accura Manual multichannel				Accura Manual XS single channel											
HTL tips	Tip volume [µl]	Cat number	0,1-2	0,5-10	1-10	2-20	5-50	10-100	20-200	100-1000	0,5-10	5-50	20-200	40-350	0,1-2	0,5-10	1-10	2-20	5-50	10-100	20-200	100-1000			
NON FILTER	0,5-10	301-03-052, 301-22-052	•	•	•	•					•				•	•	•	•							
		301-03-102, 301-22-102	•	•	•	•						•				•	•	•	•						
		301-03-152, 301-22-152	•	•	•	•						•				•	•	•	•						
		301-03-351	•	•	•	•						•				•	•	•	•						
		301-79-102	•	•	•	•						•				•	•	•	•						
	1-200	301-02-052, 301-06-052				•			•										•			•			
		301-02-152				•			•										•			•			
		301-02-202, 301-06-352				•			•										•			•			
		301-06-252				•			•										•			•			
		301-79-052				•			•										•			•			
1-300	301-16-052							•										•			•				
	301-16-102							•										•			•				
100-1000	301-01-202																						•		
	301-01-302																						•		
	301-01-352																						•		
FILTER	0,5-10	302-05-102, 302-26-152	•																						
	0,5-20	302-03-102	•			•																			
	1-100	302-08-102	•					•																	
	1-200	302-04-102																						•	
	1-300	302-12-102																							
	100-1000	302-01-102																						•	

• compatible

• incompatible

November 2020