

LumiPCR™ / ThinPCR™**0.2 mL, 8 - strip PCR tubes with individually attached flat cap**

ET080.2F5	0.2 mL, 8 - strip PCR tubes with individually attached frosted flat cap
ET080.2F35	0.2 mL, 8 - strip PCR tubes, individually attached optically clear flat caps
ET080.2F35.RW	0.2 mL, 8 - strip PCR tubes, individually attached optically clear flat caps with recessed window

Composition

	ET080.2F5	ET080.2F35, ET080.2F35.RW
Color	Clear / frosted	Optically clear caps for qPCR
Material	Medical - grade virgin polypropylene	

Storage and Usability

Storage	Protect from direct sunlight and UV light Store dry at room temperature.
Operational temperature	- 20°C to + 105°C
Usage	Designed for single use in PCR applications. Intended for research purposes only.

Endotoxin (Pyrogen) Testing: The product is tested to be in accordance with ANSI/AAMI ST 72:2002| (R)2010. Bacterial endotoxins test, Methodologies, routine monitoring and alternatives to batch testing.

RNase / DNase Testing: This product is tested and is free of any detectable RNase / DNase contamination.

Bovine Spongiform Encephalopathy (BSE) and Transmissible Spongiform Encephalopathy (TSE): None of these products or raw materials are manufactured from or come into contact with animal material.